
# BIOPSY CATALOGUE

**BFB**<sup>TM</sup>  
m e d i c a

# INDEX


<b>COMPANY PROFILE</b>	pag. 4
<b>ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE</b>	
TOTALLY REMOVE™	pag. 8
FULLY REMOVE™	pag. 10
UNLUX™ SYSTEM	pag. 12
<b>BONE MARROW BIOPSY NEEDLE - JAMSHIDI STYLE</b>	
OSTEOBELL "T"™	pag. 14
<b>BONE MARROW ASPIRATION NEEDLE</b>	
ILIAC MARROW™	pag. 15
STERNOBELL™	pag. 16
<b>BONE MARROW HARVEST NEEDLE</b>	
STERNOBELL ESPIANTO™	pag. 17
OSTEOBELL ESPIANTO™	pag. 18
<b>SOFT TISSUE BIOPSY DEVICES</b>	
ESTER™: Automatic device	pag. 19
SPEEDYBELL™: Semi-automatic device	pag. 22
SPEEDYBELL DOPPIA CORSA™: Adjustable notch semi-automatic device	pag. 24
BIOSOFT™: Automatic aspiration biopsy needle	pag. 26
CHIBELL™: Aspiration needle for soft tissue biopsy and prenatal diagnosis	pag. 28
GAMMA™: Aspiration needle for soft tissue biopsy and prenatal diagnosis	pag. 29
MOON™: Breast lesions localization needle	pag. 30

# OUR COMPANY

**BPB MEDICA™ is an Italian manufacturing company** specializing in the design, production and marketing of high-quality healthcare products for medical use and medical-surgery devices.

BPB MEDICA™ was founded in 1999 by the Bellini family, boasting thirty years of experience in the biomedical sector. The founder, Carlo Bellini Sr., started the business in 1968 and has passed down ethics, integrity and spirit of sacrifice to his heirs. Today BPB MEDICA™ has leveraged its 50 years of experience to develop new innovative product lines, growing the company on the international level.

**BPB MEDICA™'s philosophy** is to grow alongside the needs of patients, doctors and hospital staff in general. Backed by the experience acquired by the company's specialized technical personnel and thanks to newly-adopted technologies, BPB MEDICA™ has quickly managed to make a name for itself in the domestic and international markets.

**OUR NUMBERS**

**OF EXPERIENCE**

**COUNTRIES SERVED**

**CUSTOMERS**

**PROCEDURES performed with our devices**

**OUR PRODUCT LINES**

**SPINE SURGERY**

**ORTHO-BIOLOGICS**

**ASSISTED REPRODUCTION**

**BIOPSY**

**INTENSIVE CARE**

**AESTHETIC**

BPB MEDICA™ provides painstaking service to its clientele and its primary aim is product quality. Our **internal Regulatory and Quality Departments** offer full-service support to our customers in the following activities:

- Quality Systems
- Regulatory Affairs
- Technical documents preparation
- Clinical experimentation
- Vigilance
- Formation
- OEM and Private Label regulatory and quality support

Our **Quality Department** conducts rigorous tests, from the raw materials to the equipment and the finished product. This allowed the company to obtain CE, ISO 13485 and the establishment registration by FDA.

Thanks to the **internal R&D Department**, BPB MEDICA™ conducts constant research in the reference pathologies with an aim to ever-better qualifying and improving its production standards and aiding the development of new products. Functional tests are performed by **R&D Department** in concert with Quality Department to demonstrate the product performance and to make sure that new products maintain their functional characteristics even in the worst case.

BPB MEDICA™ is an established and experienced contract medical equipment manufacturer, with:

- Advanced ISO 8 cleanroom facility for manual & automatic assembly and packaging.
- Metal refinishing department for cutting, grinding, sharpening, cleaning, echogenic marking, sealing, and reduction.
- Moulding department with vertical and horizontal moulding machines.
- Computerised warehouse with top-selling items always available in stock ready to be shipped within 24 hours.

## OUR SERVICES

### QUALITY & REGULATORY DEP.

Our primary aim is product quality together with painstaking service. Our regulatory and quality teams conduct rigorous tests to this purpose.

### OEM & PRIVATE LABEL SERVICES

À la carte production, with customer's brand name and custom colour. Full-service support of Regulatory and Quality documentation.

### RESEARCH & DEVELOPMENT

Constant research to increase production and quality standards while developing new products.

### MARKETING SUPPORT

Video tutorial, case report, presence at the major medical congresses, organization of training and courses.

### PRODUCTION PROCESS CONTROL

Complete manufacturing process carried out internally, from design to final packaging.

### FOUR WEEKS DELIVERY

Thanks to the optimization of the production process, we satisfy our customer's orders within 4 weeks.

## CUTTING, GRINDING & ELECTROPOLISHING

## CLEANROOM

## MOULDING DEPARTMENT

# TOTALLY REMOVE™

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

## **BONE MARROW BIOPSY NEEDLE**

The TOTALLY REMOVE™ needle allows users to perform bone marrow biopsies on the iliac crest. With this particular system, called "TOTALLY", it is possible to avoid the dislocation manoeuvre when the sample is being withdrawn. The sample is blocked inside the system. With a single hole in the patient, it is possible to perform both a histological and a cytological sampling.

### **FEATURES:**

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

### **SURGICAL TECHNIQUE:**

1. Insert the needle, with a stylet, perpendicular to the bone surface.  
Progress forward by rotating the needle to the right and left. When the resistance decreases, the medullary cavity has been entered.
2. Remove the stylet and advance into the medullary cavity for about 2 cm.
3. Insert the TOTALLY cannula and rotate the entire device 360° twice, cutting a fragment of bone marrow tissue from the distal part, without damaging it and avoiding its dislocation.
4. Remove the TOTALLY cannula from the needle and insert a special extractor into the same TOTALLY cannula, enabling the withdrawal of the sample to be analysed.
5. If bone marrow blood needs to be drawn, attach a syringe with a Luer-lock connector to the Luer-lock connection of the needle that is still inserted in the patient.
6. At the end of the procedure, remove the needle from the patient.

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
# FULLY REMOVE™

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

## **BONE MARROW BIOPSY NEEDLE**

The FULLY REMOVE™ needle allows users to perform bone marrow biopsies on the iliac crest. With this particular system, called “FULLY”, it is possible to avoid the dislocation manoeuvre when the sample is being withdrawn.

### **FEATURES:**

- The “Trocar” tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

### **SURGICAL TECHNIQUE:**

1. Insert the needle, with a stylet, perpendicular to the bone surface.  
Progress forward by rotating the needle to the right and left. When the resistance decreases, the medullary cavity has been entered.
2. Remove the stylet and advance into the medullary cavity for about 2 cm.
3. Insert the FULLY REMOVE™ cannula and rotate the entire device 360° twice, cutting a fragment of bone marrow tissue from the distal part, without damaging it and avoiding its dislocation.
4. Extract the entire cannula needle FULLY REMOVE™ - stylet system from the patient.
5. Remove the FULLY REMOVE™ stylet from the cannula. To deposit the sample on the microscope slide, insert the finger protection device and the extractor on the opposite side from the needle handle.

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# UNLUX™ SYSTEM

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

## **ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE**

The UNLUX™ SYSTEM needle allows users to perform bone marrow biopsies on the iliac crest. With this particular system, called "UNLUX™", it is possible to avoid the dislocation manoeuvre when the sample is being withdrawn. Due to its particular features, specific technical aspects of this device have been patented.

### **FEATURES:**

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

### **SURGICAL TECHNIQUE:**

1. Insert the needle, with a stylet, perpendicular to the bone surface.  
Progress forward by rotating the needle to the right and left. When the resistance decreases, the medullary cavity has been entered.
2. Remove the stylet and advance into the medullary cavity for about 2 cm.
3. Insert the UNLUX™ cannula and rotate the entire device 360° twice, cutting a fragment of bone marrow tissue from the distal part, without damaging it and avoiding its dislocation.
4. Extract the entire cannula needle UNLUX™ - stylet system from the patient.
5. Remove the UNLUX™ stylet from the cannula. To deposit the sample on the microscope slide, insert the finger protection device and the extractor on the opposite side from the needle handle.

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## TROCAR-TIP NEEDLE FOR BONE MARROW BIOPSY

### TROCAR-TIP NEEDLE FOR BONE MARROW BIOPSY

The OSTEObELL "T"™ needle allows users to perform bone marrow biopsies on the iliac crest.

#### FEATURES:

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

The OSTEObELL "T"™ Set includes a dedicated extractor for quick ejection of the sample to be analysed, a specific guide device to protect the operator's hands during sample extraction, and a protective cap.

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## BONE MARROW ASPIRATION NEEDLE

### **BONE MARROW ASPIRATION NEEDLE**

The ILIAC MARROW™ needle is a single-use device for bone marrow aspiration from the iliac crest.

#### **FEATURES:**

- The ergonomic handle, with Luer-lock connector, and the perfectly sealed connection between the stylet and the cannula enable quick and effective bone marrow aspiration.
- The advancement system enables precise adjustment of the needle length.
- Due to the innovative "removable" advancement system, the total length of the needle can be used.
- The triple-sharpened cannula enables an easy and non-traumatic penetration.

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### **ORDER GUIDE - ILIAC MARROW™**

Gauge	Diameter (mm)	Product Code	Needle size	Unscrewed (mm)	Screwed (mm)	Only support (mm)	Only handle (mm)	Pieces for box
15G	1,80	IM1507	15G x 7cm	9	25	32	71	10

# STERNOBELL™

## BONE MARROW ASPIRATION NEEDLE

### BONE MARROW ASPIRATION NEEDLE

The STERNOBELL™ needle is a single-use device for bone marrow aspiration from the sternum and the iliac crest.

#### FEATURES:

- Ergonomic handle.
- The millimetric advancement ring facilitates the adjustment of the length of the needle, which indicates the achievable depth.
- The triple-sharpened cannula enables an easy and non-traumatic penetration.
- With the universal Luer-lock connector, any syringe can be connected.
- It is possible to completely unscrew the millimetric advancement system to make use of the total length of the needle.

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


## BONE MARROW HARVEST NEEDLE

### **BONE MARROW HARVEST NEEDLE**

STERNOBELL EXPLANT™ is a bone marrow harvest needle designed for autologous bone marrow harvest-transplantation procedures.

#### **FEATURES:**

- The two holes on the cannula are offset from each other for a quicker and greater aspiration of bone marrow and hematopoietic cells.
- Ergonomic handle.
- The millimetric advancement ring facilitates the adjustment of the length of the needle, which indicates the achievable depth.
- The triple-sharpened cannula enables an easy and non-traumatic penetration.
- With the universal Luer-lock connector, any syringe can be connected.
- It is possible to completely unscrew the millimetric advancement system to make use of the total length of the needle.

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## BONE MARROW HARVEST NEEDLE

### **BONE MARROW HARVEST NEEDLE**

OSTEOBELL EXPLANT™ is a bone marrow harvest needle designed for autologous bone marrow harvest transplantation procedures.

#### **FEATURES:**

- The five holes on the cannula are offset from each other for a quicker and greater aspiration of bone marrow and hematopoietic cells.
- The ergonomic handle improves grip and increases the precision of the biopsy procedure.
- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- With the universal Luer-lock connector, any syringe can be connected.
- The OSTEObELL EXPLANT™ Kit includes:
  - a supplementary stylet which, once inserted into the cannula, allows an additional procedure to be carried out;
  - a closing cap for protecting the user.

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### **ORDER GUIDE - OSTEObELL EXPLANT™**

Gauge	Diameter (mm)	Product Code	Needle size	Pieces for box
11G	3,00	OBE1111TC	11G x 11cm	10

# ESTER™

## AUTOMATIC DEVICE FOR SOFT TISSUE BIOPSY

### **AUTOMATIC DEVICE FOR SOFT TISSUE BIOPSY**

ESTER™ is a fully automatic, single-use device, particularly indicated for prostatic, hepatic, renal and mammary biopsies. Due to its particular features, **specific technical aspects of this device have been patented.**

#### **FEATURES:**

1. A single button that loads the cannula and the stylet.  
With this system the device can be used with one hand, ensuring greater convenience and comfort, as well as being able to use the ultrasound scan with the free hand.
2. Two possible firing buttons: one on the side and one on the back, which can be chosen by the doctor according to his/her specific manual preference.
3. The fire can take place in two distinct alternatives, chosen by the Doctor:
  1. Delay system: pressing the button lightly allows the stylet only to advance, thus exposing the sampling cavity out of the cannula. Subsequently, pressing the button until it stops allows the cannula to advance. In this way, the fire takes place through two distinct actions and the time between one and the other is defined by the doctor.
  2. Automatic system: pressing the button until it stops advances the stylet and the cannula simultaneously (the advancement time is slightly out of phase, but the fire, in this case, takes place through a single action).

ESTER™ can be equipped with a coaxial needle (Introducer) to the device, recommended in biopsies with multiple samples (egg. mammary biopsies).

ESTER™ is available in various versions:

- In the "Double Step" version: automatic and delay system
- In the "Single Step" version: automatic single shot.

#### **SURGICAL TECHNIQUE:**

1. **LOADING:** A single button that loads the cannula and the stylet.  
Pressing the button once loads the cannula, and a second time loads the stylet.
2. **INSERTION OF THE NEEDLE:** The loaded needle is inserted into the selected point for biopsy.
3. **EXPOSING OF THE SAMPLE CAVITY:** The needle is ready for the biopsy procedure: two possible firing buttons, one on the side and one at the back, and can be chosen by the doctor.  
The fire can take place in two distinct alternatives:  
Delay system: the selected button is pressed once lightly to advance the stylet and then the button is pressed until it stops to advance the cannula.  
Automatic system: pressing the button until it stops moves the cannula and the stylet forward at the same time.
4. **FIRE:** The sample is collected inside the needle. Finish by gently extracting the needle from the tissues.
5. **SAMPLE COLLECTION:** To recover the material, repeat the loading procedure indicated in point 1, and press the load button only once to uncover the sample area. Then collect the sample for the laboratory exam.

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# SPEEDYBELL™

## SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY

### **SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY**

SPEEDYBELL™ is a semi-automatic guillotine needle for soft tissue biopsy.

#### **FEATURES:**

- The echogenic marker on the end of the cannula allows constant monitoring of the needle position during a biopsy.
- The centimetre depth markings on the cannula let users control the depth reached.
- The ring handle of the device allows it to be gripped and held easily during the biopsy procedure.
- The SPEEDYBELL™ needle can be equipped with a coaxial needle (Introducer) to the device, recommended in biopsies with multiple samples (e.g. mammary biopsies).

#### **SURGICAL TECHNIQUE:**

1. Load the needle by pulling the rear knob. The needle is loaded when the locking click of the spring is heard.
2. Insert the loaded needle into the selected biopsy point, being careful not to press the rear knob.
3. When the needle is in the selected biopsy point, press and advance the rear knob manually until the final click is heard to collect the sample to be analysed.
4. Gently extract the needle from the tissue. To recover the material, repeat the loading operation by pulling the rear knob until the locking click of the spring is heard. Push the knob slowly until the sampling area is uncovered then collect the sample to be analysed.

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# SPEEDYBELL DOPPIA CORSA™

SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY WITH 10MM OR 20MM COLLECTION AREA

## **SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY WITH 10MM OR 20MM COLLECTION AREA**

SPEEDYBELL DOPPIA CORSA™ is a semi-automatic guillotine needle for soft tissue biopsy.

### **FEATURES:**

- The distinctive feature of this needle is the possibility of choosing between a 10mm or 20mm sample cavity.
- The echogenic marker on the end of the cannula allows constant monitoring of the needle position during a biopsy.
- The centimetre depth markings on the cannula let users control the depth reached.
- The ring handle of the device allows it to be gripped and held easily during the biopsy procedure.
- The SPEEDYBELL DOPPIA CORSA™ needle can be equipped with a coaxial needle (Introducer) to the device, recommended in biopsies with multiple samples (e.g. mammary biopsies).

### **SURGICAL TECHNIQUE:**

1. Load the needle by pulling the rear knob. The needle is loaded when the locking click of the spring is heard. The first click indicates a 10mm sampling area, the second click indicates a 20mm sampling area.
2. Insert the loaded needle into the selected biopsy point, being careful not to press the rear knob.
3. When the needle is in the selected biopsy point, press and advance the rear knob manually until the final click is heard to collect the sample to be analysed.
4. Gently extract the needle from the tissue. To recover the material, repeat the loading operation by pulling the rear knob until the locking click of the spring is heard. Push the knob slowly until the sampling area is uncovered then collect the sample to be analysed.

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
# BIOSOFT™

## AUTOMATIC ASPIRATION DEVICE FOR HISTOLOGICAL BIOPSY

### **AUTOMATIC ASPIRATION DEVICE FOR HISTOLOGICAL BIOPSY**

BIOSOFT™ is a disposable automatic aspiration device used in the main histological biopsies.

#### **FEATURES:**

- Possibility to perform the biopsy with one hand, thanks to the presence of a spring that allows the sample to be taken automatically.
- Speed and certainty of biopsy performance.
- The particular sharpening of the tip enables users to obtain an integral sample, through a simple and non-traumatic penetration.
- The centimetre depth markings allow a safe positioning of the needle.
- The sliding stopper enables pre-definition of the penetration depth.
- The internal echogenic marker prevents the dragging of tumour cells during the biopsy.
- Two versions available:
  - With a fixed cannula
  - With a removable cannula

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# CHIBELL™

## ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

### ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

CHIBELL™ can be used for various types of biopsies and prenatal diagnosis: cytological aspiration, amniocentesis, chorionic villus sampling, cystocentesis, and needle aspirations.

#### FEATURES:

- The particular sharpening of the Chiba-tip improves penetration, making it less invasive and traumatic.
- The centimetre depth marker, sliding stopper and internal echo marker make it easier and more accurate to control the depth reached.
- Equipped with a Luer-lock universal connection.
- The passage of liquids can be seen through the transparent plastic handle.

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# GAMMA™

## ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

### **ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS**

GAMMA™ can be used for various types of biopsies and prenatal diagnosis: cytological aspiration, amniocentesis, chorionic villus sampling, cystocentesis, and needle aspirations.

#### **FEATURES:**

- The particular sharpening of the modified Chiba-tip improves penetration, making it less invasive and traumatic.
- The centimetre depth marker, sliding stopper and internal echo marker make it easier and more accurate to control the depth reached.
- Equipped with a Luer-lock universal connection.
- The passage of liquids can be seen through the transparent plastic handle.

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# MOON™

## NEEDLE FOR THE LOCALIZATION OF BREAST LESIONS

### **NEEDLE FOR THE LOCALIZATION OF BREAST LESIONS**

MOON™ is an effective support in the diagnostic imaging technique for the preoperative localization of non-palpable breast lesions.

#### **FEATURES:**


- The particular sharpness of the needle tip improves penetration.
- The centimetre depth marker, sliding stopper and internal echo marker, make it easier and more accurate to control the depth reached.
- The non-repositionable hook, available with a “single hook” or “double hook”, is used as a guide towards the lesion to be eliminated.

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FDA (Establishment Registration Number) 9617616

©Copyright 2022 - All Rights Reserved

Rev. 03 del 31-05-2022

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Società soggetta ad attività di direzione e coordinamento da parte della società Bpunto3 S.r.l.