
BIOPSY CATALOGUE

3PBTM
medica

INDEX

COMPANY PROFILE	pag. 3
ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE	
TOTALLY REMOVE	pag. 4
FULLY REMOVE	pag. 6
UNLUX SYSTEM	pag. 8
BONE MARROW BIOPSY NEEDLE - JAMSHIDI STYLE	
OSTEOBELL "T"	pag. 10
BONE MARROW ASPIRATION NEEDLE	
ILIAC MARROW	pag. 11
STERNOBELL	pag. 12
BONE MARROW HARVEST NEEDLE	
STERNOBELL ESPIANTO	pag. 13
OSTEOBELL ESPIANTO	pag. 14
SOFT TISSUE BIOPSY DEVICES	
ESTER: Automatic device	pag. 15
SPEEDYBELL: Semi-automatic device	pag. 18
SPEEDYBELL DOPPIA CORSA: Adjustable notch semi-automatic device	pag. 20
BIOSOFT: Automatic aspiration biopsy needle	pag. 22
CHIBELL: Aspiration needle for soft tissue biopsy and prenatal diagnosis	pag. 24
GAMMA: Aspiration needle for soft tissue biopsy and prenatal diagnosis	pag. 25
MOON: Breast lesions localization needle	pag. 26

The Bellini family has been operating in the biomedical sector for three generations, for over 50 years. The founder, Carlo Bellini, started the business in 1968 and passed ethics, integrity and a spirit of sacrifice down to his heirs. Today BPB MEDICA™ continues to manufacture all its components and devices within the company, with sales being made through our collaborator-salespeople in Italy and the rest of the world.

The philosophy of BPB MEDICA™ is to grow alongside the needs of patients, doctors and hospital staff in general. BPB MEDICA™ was founded in 1999 by partners boasting thirty year's experience in the biomedical sector. The company is engaged in the design, production and marketing of healthcare products for medical use and medical-surgical devices. Thanks to both the experience gained by the specialized technical staff working in-house as well as the adoption of new technologies, BPB MEDICA™ has managed to establish itself rapidly on the global market.

BPB MEDICA™ operates with state-of-the-art production machinery and equipment and the entire production process is carried out in-house (from the design to the final packaging).

BPB MEDICA™ provides painstaking service to its clientele and its primary aim is product quality. The rigorous tests carried out for this purpose, from the raw materials to the equipment and the finished product, have enabled the company to obtain CE, ISO 13485 and the establishment registration by the FDA.

BPB MEDICA™ carries out constant research and collaborates with highly specialized professionals in the sector in order to increase production and quality standards and develop new products.

TOTALLY REMOVE

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

BONE MARROW BIOPSY NEEDLE

The TOTALLY REMOVE needle allows users to perform bone marrow biopsies on the iliac crest. With this particular system, called "TOTALLY", it is possible to avoid the dislocation manoeuvre when the sample is being withdrawn. The sample is blocked inside the system. With a single hole in the patient, it is possible to perform both a histological and a cytological sampling.

FEATURES:

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-Lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

SURGICAL TECHNIQUE:

1. Insert the needle, with a stylet, perpendicular to the bone surface.
Progress forward by rotating the needle to the right and left. When the resistance decreases, the medullary cavity has been entered.
2. Remove the stylet and advance into the medullary cavity for about 2 cm.
3. Insert the TOTALLY cannula and rotate the entire device 360° twice, cutting a fragment of bone marrow tissue from the distal part, without damaging it and avoiding its dislocation.
4. Remove the TOTALLY cannula from the needle and insert a special extractor into the same TOTALLY cannula, enabling the withdrawal of the sample to be analysed.
5. If bone marrow blood needs to be drawn, attach a syringe with a Luer-Lock connector to the Luer-Lock connection of the needle that is still inserted in the patient.
6. At the end of the procedure, remove the needle from the patient.

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FULLY REMOVE

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

BONE MARROW BIOPSY NEEDLE

The FULLY REMOVE needle allows users to perform bone marrow biopsies on the iliac crest. With this particular system, called "FULLY", it is possible to avoid the dislocation manoeuvre when the sample is being withdrawn.

FEATURES:

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-Lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

SURGICAL TECHNIQUE:

1. Insert the needle, with a stylet, perpendicular to the bone surface.
Progress forward by rotating the needle to the right and left. When the resistance decreases, the medullary cavity has been entered.
2. Remove the stylet and advance into the medullary cavity for about 2 cm.
3. Insert the FULLY cannula and rotate the entire device 360° twice, cutting a fragment of bone marrow tissue from the distal part, without damaging it and avoiding its dislocation.
4. Extract the entire cannula needle FULLY - stylet system from the patient.
5. Remove the FULLY stylet from the cannula. To deposit the sample on the microscope slide, insert the finger protection device and the extractor in the opposite side from the needle handle.

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UNLUX SYSTEM

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

ANTI-DISLOCATION BONE MARROW BIOPSY NEEDLE

The UNLUX SYSTEM needle allows users to perform bone marrow biopsies on the iliac crest. With this particular system, called "UNLUX", it is possible to avoid the dislocation manoeuvre when the sample is being withdrawn.

Due to its particular features, specific technical aspects of this device have been patented.

FEATURES:

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-Lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

SURGICAL TECHNIQUE:

1. Insert the needle, with a stylet, perpendicular to the bone surface.
Progress forward by rotating the needle to the right and left. When the resistance decreases, the medullary cavity has been entered.
2. Remove the stylet and advance into the medullary cavity for about 2 cm.
3. Insert the UNLUX cannula and rotate the entire device 360° twice, cutting a fragment of bone marrow tissue from the distal part, without damaging it and avoiding its dislocation.
4. Extract the entire cannula needle UNLUX - stylet system from the patient.
5. Remove the UNLUX stylet from the cannula. To deposit the sample on the microscope slide, insert the finger protection device and the extractor in the opposite side from the needle handle.

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TROCAR-TIP NEEDLE FOR BONE MARROW BIOPSY

TROCAR-TIP NEEDLE FOR BONE MARROW BIOPSY

The OSTEOBELL "T" needle allows users to perform bone marrow biopsies on the iliac crest.

FEATURES:

- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- The ergonomic handle improves grip and the precision of the biopsy operation.
- With the universal Luer-Lock connector, any syringe can be connected.
- The ultra-sharp crown tip of the cannula facilitates the removal of an intact bone marrow fragment.

The "Osteobell T" Set includes: a dedicated extractor for quick ejection of the sample to be analysed; a specific guide device to protect the operator's hands during sample extraction; a protective cap.

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BONE MARROW ASPIRATION NEEDLE

BONE MARROW ASPIRATION NEEDLE

The ILIAC MARROW needle is a single-use device for bone marrow aspiration from the iliac crest.

FEATURES:

- The ergonomic handle, with Luer-Lock connector, and the perfectly sealed connection between the stylet and the cannula enables quick and effective bone marrow aspiration.
- The advancement system enables precise adjustment of the needle length.
- Due to the innovative "removable" advancement system, the total length of the needle can be used.
- The triple sharpened cannula enables an easy and non-traumatic penetration.

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ORDER GUIDE - ILIAC MARROW

Gauge	Diameter (mm)	Product Code	Needle size	Unscrewed (mm)	Screwed (mm)	Only support (mm)	Only handle (mm)	Pieces for box
15G	1,80	IM1507	15G x 7cm	9	25	32	71	10

STERNOBELL

BONE MARROW ASPIRATION NEEDLE

BONE MARROW ASPIRATION NEEDLE

The STERNOBELL needle is a single-use device for bone marrow aspiration from the sternum and the iliac crest.

FEATURES:

- Ergonomic handle.
- The millimetric advancement ring facilitates the adjustment of the length of the needle, which indicates the achievable depth.
- The triple sharpened cannula enables an easy and non-traumatic penetration.
- With the universal Luer-Lock connector, any syringe can be connected.
- It is possible to completely unscrew the millimetric advancement system to make use of the total length of the needle.

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
BONE MARROW HARVEST NEEDLE

BONE MARROW HARVEST NEEDLE

STERNOBELL ESPIANTO is a bone marrow harvest needle designed for allogeneic and autologous bone marrow harvest-transplantation procedures.

FEATURES:

- The two holes on the cannula are offset from each other for a quicker and greater aspiration of bone marrow and hematopoietic cells.
- Ergonomic handle.
- The millimetric advancement ring facilitates the adjustment of the length of the needle, which indicates the achievable depth.
- The triple sharpened cannula enables an easy and non-traumatic penetration.
- With the universal Luer-Lock connector, any syringe can be connected.
- It is possible to completely unscrew the millimetric advancement system to make use of the total length of the needle.

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BONE MARROW HARVEST NEEDLE

BONE MARROW HARVEST NEEDLE

OSTEOBELL ESPIANTO is a bone marrow harvest needle designed for allogeneic and autologous bone marrow transplantations procedures.

FEATURES:

- The five holes on the cannula are offset from each other for a quicker and greater aspiration of bone marrow and hematopoietic cells.
- The ergonomic handle improves grip, and increases the precision of the biopsy procedure.
- The "Trocar" tip stylet guarantees fast and non-traumatic penetration.
- With the universal Luer-Lock connector, any syringe can be connected.
- The OSTEOBELL ESPIANTO Kit includes:
 - an supplementary stylet which, once inserted into the cannula, allows an additional procedure to be carried out;
 - a closing cap for protecting the user.

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ORDER GUIDE - OSTEOBELL ESPIANTO

Gauge	Diameter (mm)	Product Code	Needle size	Pieces for box
11G	3,00	OBE1111TC	11G x 11cm	10

ESTER

AUTOMATIC DEVICE FOR SOFT TISSUE BIOPSY

AUTOMATIC DEVICE FOR SOFT TISSUE BIOPSY

ESTER is a fully automatic, single-use device, particularly indicated for prostatic, hepatic, renal and mammary biopsies. Due to its particular features, **specific technical aspects of this device have been patented**.

FEATURES:

1. A single button that loads the cannula and the stylet.

With this system the device can be used with one hand, ensuring greater convenience and comfort, as well being able to use of the ultrasound scan with the free hand.

2. Two possible firing buttons: one on the side and one on the back, which can be chosen by the doctor in relation to his/her specific manual preference.

3. The fire can take place in two distinct alternatives, chosen by the Doctor:

3.1 Delay system: pressing the button lightly allows the stylet only to advance, thus exposing the sampling cavity out of the cannula.

Subsequently, pressing the button until it stops allows the cannula to advance forward. In this way the fire takes place through two distinct actions and the time between one and the other is defined by the doctor.

3.2 Automatic system: pressing the button until it stops advances the stylet and the cannula simultaneously (the advancement time is slightly out of phase, but the fire, in this case, takes place through a single action).

ESTER can be equipped with a coaxial needle (Introducer) to the device, recommended in biopsies with multiple samples (egg. mammabiopsies).

ESTER is available in various versions:

- In the "Double Step" version: automatic and delay system
- In the "Single Step" version: automatic single shot.

SURGICAL TECHNIQUE:

1. LOADING: A single button that loads the cannula and the stylet.

Pressing the button once loads the cannula, and a second time loads the stylet.

2. INSERTION OF THE NEEDLE: The loaded needle is inserted in the selected point for biopsy.

3. EXPOSING OF THE SAMPLE CAVITY: The needle is ready for the biopsy procedure: two possible firing buttons, one on the side and one at the back, and can be chosen by the doctor.

The fire can take place in two distinct alternatives:

Delay system: the selected button is pressed once lightly to advance the stylet and then the button is pressed until it stops to advance the cannula.

Automatic system: pressing the button until it stops moves the cannula and the stylet forward at the same time.

4. FIRE: The sample is collected inside the needle. Finish by gently extracting the needle from the tissues.

5. SAMPLE COLLECTION: To recover the material, repeat the loading procedure indicated in point 1, press the load button only once to uncover the sample area. Then collect the sample for the laboratory exam.

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SPEEDYBELL

SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY

SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY

SPEEDYBELL is a semi-automatic guillotine needle for soft tissue biopsy.

FEATURES:

- The echogenic marker on the end of the cannula allows constant monitoring of the needle position during biopsy.
- The centimetre depth markings on the cannula let users control of the depth reached.
- The ring handle of the device allows it to be gripped and held easily during the biopsy procedure.
- The SPEEDYBELL needle can be equipped with a coaxial needle (Introducer) to the device, recommended in biopsies with multiple samples (e.g. mammary biopsies).

SURGICAL TECHNIQUE:

1. Load the needle by pulling the rear knob. The needle is loaded when the locking click of the spring is heard.
2. Insert the loaded needle into the selected biopsy point, being careful not to press the rear knob.
3. When the needle is in the selected biopsy point, press and advance the rear knob manually until the final click is heard to collect the sample to be analysed.
4. Gently extract the needle from the tissue. To recover the material, repeat the loading operation by pulling the rear knob until the locking click of the spring is heard. Push the knob slowly until the sampling area is uncovered then collect the sample to be analysed.

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SPEEDYBELL DOPPIA CORSA

SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE
BIOPSY WITH 10MM OR 20MM COLLECTION AREA

SEMI-AUTOMATIC GUILLOTINE NEEDLE FOR SOFT TISSUE BIOPSY WITH 10MM OR 20MM COLLECTION AREA

SPEEDYBELL DOPPIA CORSA is a semi-automatic guillotine needle for soft tissue biopsy.

FEATURES:

- The distinctive feature of this needle is the possibility of choosing between a 10mm or 20mm sample cavity.
- The echogenic marker on the end of the cannula allows constant monitoring of the needle position during biopsy.
- The centimetre depth markings on the cannula let users control of the depth reached.
- The ring handle of the device allows it to be gripped and held easily during the biopsy procedure.
- The SPEEDYBELL DOPPIA CORSA needle can be equipped with a coaxial needle (Introducer) to the device, recommended in biopsies with multiple samples (e.g. mammary biopsies).

SURGICAL TECHNIQUE:

1. Load the needle by pulling the rear knob. The needle is loaded when the locking click of the spring is heard. The first click indicates a 10mm sampling area, the second click indicates a 20mm sampling area.
2. Insert the loaded needle into the selected biopsy point, being careful not to press the rear knob.
3. When the needle is in the selected biopsy point, press and advance the rear knob manually until the final click is heard to collect the sample to be analysed.
4. Gently extract the needle from the tissue. To recover the material, repeat the loading operation by pulling the rear knob until the locking click of the spring is heard. Push the knob slowly until the sampling area is uncovered then collect the sample to be analysed.

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BIOSOFT

AUTOMATIC ASPIRATION DEVICE FOR HISTOLOGICAL BIOPSY

AUTOMATIC ASPIRATION DEVICE FOR HISTOLOGICAL BIOPSY

BIOSOFT is a disposable automatic aspiration device used in the main histological biopsies.

FEATURES:

- Possibility to perform the biopsy with one hand, thanks to the presence of a spring that allows the sample to be taken automatically.
- Speed and certainty of biopsy performance.
- The particular sharpening of the tip enables users to obtain an integral sample, through a simple and non-traumatic penetration.
- The centimetre depth markings allow a safe positioning of the needle.
- The sliding stopper enables pre-definition of the penetration depth.
- The internal echogenic marker prevents dragging of tumour cells during biopsy.
- Two versions available:
 - With a fixed cannula
 - With a removable cannula

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CHIBELL

ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

CHIBELL can be used for various types of biopsies and prenatal diagnosis: cytological aspiration, amniocentesis, chorionic villus sampling, cystocentesis, and needle aspirations.

FEATURES:

- The particular sharpening of the Chiba-tip improves penetration, making it less invasive and traumatic.
- The centimetre depth marker, sliding stopper and internal echo marker, make it easier and more accurate to control the depth reached.
- Equipped with a Luer-Lock universal connection.
- The passage of liquids can be seen through the transparent plastic handle.

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GAMMA

ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

ASPIRATION NEEDLE FOR SOFT TISSUE BIOPSY AND PRENATAL DIAGNOSIS

GAMMA can be used for various types of biopsies and prenatal diagnosis: cytological aspiration, amniocentesis, chorionic villus sampling, cystocentesis, and needle aspirations.

FEATURES:

- The particular sharpening of the modified Chiba-tip improves penetration, making it less invasive and traumatic.
- The centimetre depth marker, sliding stopper and internal echo marker, make it easier and more accurate to control the depth reached.
- Equipped with a Luer-Lock universal connection.
- The passage of liquids can be seen through the transparent plastic handle.

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MOON

NEEDLE FOR THE LOCALIZATION OF BREAST LESIONS

NEEDLE FOR THE LOCALIZATION OF BREAST LESIONS

MOON is an effective support in the diagnostic imaging technique for the preoperative localization of non-palpable breast lesions.

FEATURES:

- The particular sharpness of the needle tip improves penetration.
- The centimetre depth marker, sliding stopper and internal echo marker, make it easier and more accurate to control the depth reached.
- The non-repositionable hook, available with a "single hook" or "double hook", is used as a guide towards the lesion to be eliminated.

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CUTTING, GRINDING & ELECTROPOOLISHING

CLEAN ROOM

MOLDING DEPARTMENT

BIOPSYBELL S.R.L.
Via Aldo Manuzio 24
41037 Mirandola (MO) - Italy
Tel. +39 0535 27850
Fax. +39 0535 33526
international1@biopsybell.it
www.biopsybell.it

FDA (Establishment Registration Number) 9617616
©Copyright 2019 - All Rights Reserved

Rev. 01 del 31-10-19