


# ASSISTED REPRODUCTION TECHNOLOGY



**BPB**<sup>™</sup>  
m e d i c a



# INDEX

<b>COMPANY PROFILE</b>	pag. 3
<b>OVUM ASPIRATION NEEDLE</b>	
KINDER PLUS™: Single Lumen	pag. 6
KINDER DL™: Double Lumen	pag. 7
<b>VITRIFICATION CARRIER</b>	
SUPERVITRI™	pag. 8

# OUR COMPANY




**BPB MEDICA™ is an Italian manufacturing company** specializing in the design, production and marketing of high-quality healthcare products for medical use and medical-surgery devices.

BPB MEDICA™ was founded in 1999 by the Bellini family, boasting thirty years of experience in the biomedical sector. The founder, Carlo Bellini Sr., started the business in 1968 and has passed down ethics, integrity and spirit of sacrifice to his heirs. Today BPB MEDICA™ has leveraged its 50 years of experience to develop new innovative product lines, growing the company on the international level.



OF EXPERIENCE



# OUR COMPANY

**BPB MEDICA's™ philosophy** is to grow alongside the needs of patients, doctors and hospital staff in general. Backed by the experience acquired by the company's specialized technical personnel and thanks to newly-adopted technologies, BPB MEDICA™ has quickly managed to make a name for itself in the domestic and international markets.



**COUNTRIES SERVED**



BPB MEDICA™ provides painstaking service to its clientele and its primary aim is product quality. The **internal Regulatory and Quality Departments** conduct rigorous tests, from the raw materials to the equipment and the finished product. This allowed the company to obtain CE, ISO 13485 and the establishment registration by FDA.



BPB MEDICA™ operates with high qualitative standards aimed at increasing customer satisfaction through a continuous improvement regime. We conduct multiple quality tests during the entire production process: this allows us to reduce the production waste and grant high-quality products.



Thanks to the **internal R&D Department** BPB MEDICA™ conducts constant research on the reference pathologies with an aim to ever-better qualifying and improve its production standards and aid the development of new products.



## OUR PRODUCT LINES



**SPINE  
SURGERY**



**ORTHO-  
BIOLOGICS**



**ASSISTED  
REPRODUCTION**



**BIOPSY**



**INTENSIVE  
CARE**



**AESTHETIC**



**CLEANROOM**  
Manual & automatic assembly, packaging.



**METAL REFINISHING DEPARTMENT**  
Cutting, grinding, sharpening, cleaning, echogenic marking, sealing, reduction.



**MOULDING DEPARTMENT**  
Vertical and horizontal moulding.

## OUR SERVICES



### **QUALITY & REGULATORY DEP.**

Our primary aim is product quality and to provide painstaking service. Our regulatory and quality teams conduct rigorous tests for this purpose.



### **OEM & PRIVATE LABEL SERVICES**

À la carte production, with customer's brand name and custom colour. Full-service support of Regulatory and Quality documentation.



### **RESEARCH & DEVELOPMENT**

Constant research to increase production and quality standards while developing new products.



### **MARKETING SUPPORT**

Video tutorial, case report, presence at the major medical congresses, organization of training and courses.



### **PRODUCTION PROCESS CONTROL**

Complete manufacturing process carried out internally, from design to final packaging.




### **PERMANENT STOCK**

Top selling items are always available in stock and can be shipped within 24 hours.



# KINDER PLUS™

## SINGLE LUMEN OVUM ASPIRATION NEEDLE



Direct connection between needle and aspiration line

Ergonomic grip with thumb notch to indicate tip direction

Echo markings for accurate placement under ultrasound guidance

Triple cut tip to reduce trauma


**DIAMETERS AND LENGTHS AVAILABLE:**

17G - 18G  
30 cm - 33 cm - 35 cm

**Different measures are available upon request**

# KINDER DL™

## DOUBLE LUMEN OVUM ASPIRATION NEEDLE



**DIAMETERS AND LENGTHS AVAILABLE:**

16G - 17G  
30 cm - 33 cm - 35 cm

**Different measures are available upon request**

# SUPERVITRI™

## VITRIFICATION CARRIER

Hermetic closure  
protects  
oocyte/embryo  
from contamination

**PATENTED**

**LOCKING  
SYSTEM**

Direction indicator for easier  
and safer cell loading

Rough surface area for  
writings and recording  
patient information


**"U" TIP**

>Thin strip to obtain a  
superfast cooling rate;  
>Curved design ("U" tip) for  
oocyte/embryo protection  
during loading and sealing;

**8**

**MAXIMUM STORAGE:  
8 SUPERVITRI™ IN VISIOTUBE**

**5 COLORS AVAILABLE**







**BIOPSYBELL S.R.L. Società Unipersonale**

Via Aldo Manuzio 24

41037 Mirandola (MO) - Italy

Tel. +39 0535 27850

Fax. +39 0535 33526

[international1@biopsybell.it](mailto:international1@biopsybell.it)

[www.biopsybell.com](http://www.biopsybell.com)

Società soggetta ad attività di direzione e coordinamento da parte della società Bpunto3 S.r.l.

FDA (Establishment Registration Number) 9617616

©Copyright 2022 - All Rights Reserved



Rev. 03 del 25-05-2022