

BEAUTY-STEM™

SUPERIOR ADIPOSE TISSUE

VOLUME RESTORATION FOR AESTHETIC & PLASTIC SURGERY

BEAUTY-STEM™ is a closed-circuit single-use kit for adipose tissue purification without any centrifuge and with minimal manipulation that allows superior and long-lasting engraftment.

The final purified product is ideal for remodelling and volume restoration.



KEY ADVANTAGES



TIME EFFICIENT

- Processes up to 400 ml of adipose tissue in 10 minutes.
- Compared to centrifugation, **BEAUTY-STEM™** allows for saving time and related costs.
- All-in-one system that processes, purifies and extracts lipoaspirate within one device.
- More than one procedure with the same kit in one surgical session.



MINIMAL MANIPULATION

- **BEAUTY-STEM™** preserves the biological properties of the cells and maximizes cells engraftment.
- Continuous saline solution washing eliminates any traumatic action that may damage the extracellular matrix.
- Sterile closed-circuit environment that zeroes the risk of air and bacteria contamination.




EASY TO USE

- No centrifugation is required.
- Minimal training, preparation and clean-up.
- Only 1 operator required.
- Ready to use autologous lipofilling.

BEAUTY-STEM™


SUPERIOR ADIPOSE TISSUE




HOW IT WORKS:

NO CENTRIFUGE. The entire processing phase of the liposuctioned tissue occurs inside the device thanks to continuous saline solution washing:

- No traumatic action that may damage the extracellular matrix and its essential trophic and anti-inflammatory function;
- Reduced cellular stress.



The collection and processing bag purifies the drawn tissue preserving its biological properties and maximizing its regenerative potential.



The end product is a compound and purified autologous adipose tissue that keeps the biological properties of the original tissue intact.

The whole extracellular matrix acts as a natural structure for the cells, increasing their vitality and assuring volume restoration.



The results of lipofilling last much longer than any other filler on the market.

APPLICATIONS



FACE

Restoration of facial volumes and youthfulness



BURNS

Improvement of the natural healing process of burned tissue



GENITAL AREA

Female auto augmentation of labia majora.
Male auto augmentation of penis girth.



SCARS

Treatment of hypertrophic scars



BODY

Remodelling and volume restoration of breast and gluteus area for aesthetic or reconstructive reasons.

