

INTRAOX™

AGO PER INFUSIONE INTRAOSSEA

È un dispositivo monouso utilizzato per la somministrazione e l'infusione intraossea di liquidi e farmaci in situazioni di emergenza, laddove non è possibile un accesso vascolare tradizionale o quando non è possibile ottenerlo in tempi brevi. L'accesso intraosseo può essere una soluzione valida e rapida nel trattamento di pazienti traumatizzati in condizioni critiche.

- I fori sulla cannula consentono una più rapida infusione ed una distribuzione multi-direzionale dei liquidi iniettati.
- L'impugnatura ergonomica migliora la presa ed aumenta la precisione dell'accesso intraosseo.
- L'attacco luer-lock universale permette la connessione con qualsiasi siringa.

INTRAOSSSEOUS INFUSION NEEDLE

It is a disposable device used to provide and infuse liquids and drugs under emergency conditions, when the common vascular access is not available or not immediately available. Intraosseous access can be an appropriate and quick alternative to treat shocked patients in critical conditions.

- *The holes on the cannula allow a more rapid infusion and a multi-directional distribution of the injected liquids.*
- *The ergonomic handle improves grip and the precision of the intraosseous access.*
- *A universal luer lock connector suitable with any kind of syringe.*



Guida all'ordine Order guide

Calibro Gauge	Diametro Diameter (mm)	3 cm	Q.tà per scatola Q.ty for box
14G	2,00	IOF140305	10
15G	1,80	IOF150305	10
16G	1,60	IOF160305	10
18G	1,20	IOF180305	10